

Percentage of Water in Popcorn

OBJECTIVES

- *Measure the masses of various combinations of a beaker, oil, and popcorn kernels.*
- *Determine the percentages of water in popcorn kernels.*
- *Compare experimental data.*

MATERIALS

- aluminum foil (1 sheet)
- beaker, 250 mL
- Bunsen burner with gas tubing and striker
- kernels of popcorn for each of three brands (80)
- oil to coat the bottom of the beaker
- ring stand, iron ring, and wire gauze

BACKGROUND

Popcorn pops because of the natural moisture inside each kernel. When the internal water is heated above 100°C, the liquid water changes to a gas, which takes up much more space than the liquid, so the kernel expands rapidly.

The percentage of water in popcorn can be determined by the following equation.

$$\frac{\text{initial mass} - \text{final mass}}{\text{initial mass}} \times 100 = \text{percentage of H}_2\text{O in unpopped popcorn}$$

The popping process works best when the kernels are first coated with a small amount of vegetable oil. Make sure you account for the presence of this oil when measuring masses. In this lab, you will design a procedure for determining the percentage of water in three samples of popcorn. The popcorn is for testing only, and *must not* be eaten.

SAFETY



For review of safety, please see **Safety in the Chemistry Laboratory** in the front of your book.

PREPARATION

1. In your notebook, prepare a data table like the one your teacher has made on the board. The table should have three columns, one for each brand of popcorn.