



# Reduction of Manganese in Permanganate Ion

## OBJECTIVES

- *Demonstrate* proficiency in performing redox titrations and recognizing end points of a redox reaction.
- *Write* a balanced oxidation-reduction equation for a redox reaction.
- *Determine* the concentration of a solution by using stoichiometry and volume data from a titration.

## MATERIALS

- 0.0200 M  $\text{KMnO}_4$
- 1.0 M  $\text{H}_2\text{SO}_4$
- 100 mL graduated cylinder
- 125 mL Erlenmeyer flasks, 4
- 250 mL beakers, 2
- 400 mL beaker
- burets, 2
- distilled water
- double buret clamp
- $\text{FeSO}_4$  solution
- ring stand
- wash bottle

## BACKGROUND

In Chapter 15, you studied acid-base titrations in which an unknown amount of acid is titrated with a carefully measured amount of base. In this procedure, a similar approach called a *redox titration* is used. In a redox titration, the reducing agent,  $\text{Fe}^{2+}$ , is oxidized to  $\text{Fe}^{3+}$  by the oxidizing agent,  $\text{MnO}_4^-$ . When this process occurs, the Mn in  $\text{MnO}_4^-$  changes from a +7 to a +2 oxidation state and has a noticeably different color. You can use this color change to signify a redox reaction “end point.” When the reaction is complete, any excess  $\text{MnO}_4^-$  added to the reaction mixture will give the solution a pink or purple color. The volume data from the titration, the known molarity of the  $\text{KMnO}_4$  solution, and the mole ratio from the balanced redox equation will give you the information you need to calculate the molarity of the  $\text{FeSO}_4$  solution.

## SAFETY



For review of safety, please see **Safety in the Chemistry Laboratory** in the front of your book.

## PREPARATION

1. In your lab notebook, prepare a data table like the one shown on the next page.
2. Clean two 50 mL burets with a buret brush and distilled water. Rinse each buret at least three times with distilled water to remove contaminants.
3. Label one 250 mL beaker “0.0200 M  $\text{KMnO}_4$ ” and the other “ $\text{FeSO}_4$ .” Label three of the flasks “1,” “2,” and “3.” Label the 400 mL beaker “Waste.” Label one buret “ $\text{KMnO}_4$ ” and the other “ $\text{FeSO}_4$ .”