

# Determining pH and Titrations

## SECTION 2


### OBJECTIVES

- Describe how an acid-base indicator functions.
- Explain how to carry out an acid-base titration.
- Calculate the molarity of a solution from titration data.

## Indicators and pH Meters

An approximate value for the pH of a solution can be obtained using acid-base indicators. **Acid-base indicators** are compounds whose colors are sensitive to pH. In other words, the color of an indicator changes as the pH of a solution changes.


Indicators change colors because they are either weak acids or weak bases. In solution, the equilibrium of an indicator that is a weak acid can be represented by the equation below, which is modeled in **Figure 4**.



( $\text{In}^-$  is the symbol of the anion part of the indicator.) The colors that an indicator displays result from the fact that  $\text{HIn}$  and  $\text{In}^-$  are different colors.

In acidic solutions, any  $\text{In}^-$  ions that are present act as Brønsted bases and accept protons from the acid. The indicator is then present in largely nonionized form,  $\text{HIn}$ . The indicator has its acid-indicating color, as shown for litmus in **Figure 4**.

In basic solutions, the  $\text{OH}^-$  ions from the base combine with the  $\text{H}^+$  ions produced by the indicator. The indicator molecules further ionize to offset the loss of  $\text{H}^+$  ions. The indicator is thus present largely in the form of its anion,  $\text{In}^-$ . The solution now displays the base-indicating color, which for litmus is blue.



**FIGURE 4** Basic solutions shift the equilibrium of litmus to the right. The ionized form,  $\text{In}^-$ , then predominates, and the litmus turns blue. Acidic solutions shift the equilibrium of the indicator litmus to the left. The nonionized form,  $\text{HIn}$ , predominates, and the litmus turns red.